


SINEAD

(Software Integrating NEuroimaging And other Data)



Wong-Lin, Kongfatt <k.wong-lin@ulster.ac.uk>, Intelligent Systems Research Centre, Ulster University
 Vahab Youssefzadeh <Vahab.Youssefzadeh@cchmc.org>, Cincinnati Children's Hospital Medical Center

June, 2016

SINEAD

This Matlab based software provides a GUI framework for conducting machine learning on imaging data and in particular for an AD application.

It integrates,

- SPM Toolbox (www.fil.ion.ucl.ac.uk/spm/) for preprocessing of imaging data,
- PRoNTO (www.mlnl.cs.ucl.ac.uk/pronto) for machine learning and combined analysis of imaging data,
- Surfstat (<http://www.math.mcgill.ca/keith/surfstat/>) and freesurfer (<http://freesurfer.net/>) for cortical thickness analysis.


References regarding SINEAD are listed below:

1. Youssofzadeh V, McGuinness B, Maguire L and Wong-Lin K, 2016. A kernel-based Gaussian process for combined MRI and PET in Alzheimer's disease. In: Alzheimer's Research UK (ARUK) 2016 Conference, Manchester, UK. 8-9 March 2016.
2. Youssofzadeh V, McGuinness B, Maguire L and Wong-Lin K. Multi-kernel learning with Dartel enhances MRI-PET classification and prediction of Alzheimer's disease: group and individual analyses, In preparation.

Installing SINEAD


1. Download CADSuM on your computer and add the directory path with subdirectories to your Matlab path.
2. Ensure you have the SPM12, PRoNTO directory in your Matlab path as well.

Pre-processing of scans





- SPM12 toolbox has to be added to Matlab routines

Set origin




Set origin to AC area was completed!




Preprocessing

- Data/scans are selected using batch interface. A sample batch is available in “SPM batch” folder.


MRI preprocessing




PET preprocessing




Preprocessing, checking results



MRI





PET



Tissue volumes

- For more details about TV, please see, Ashburner 2015 VBM Tutorial that is provided in the Manual folder!




Volume 1: gray matter

Volume 2: white matter

Volume 3: CSF

Total intracranial volume = Volume 1 + Volume 2 + Volume 3




Machine learning



- PRoNTo toolbox is utilized for the machine learning process.


Data Selection

- Data/scans are selected using batch interface.



A sample batch is available in "PRoNTTo batch" folder. It is related to 3*20 scans from 3 groups of AD, MCI and Healthy elderly subjects .


Building kernels



3




```
SPM12: spm_reslice (v5929) 11:16:03 - 26/06/2016
=====
Completed : 11:16:03 - 26/06/2016
Feature set 'Fset' not found in PRT.mat. Creating...
File array does not exist for modality 'GM'. Creating...
File array does not exist for modality 'PiB'. Creating...
Saving feature set to: PRT.mat.....>>
Saving kernel to: Fset.mat.....>>
Done.
fx >> |
```

An example of combined PET–MRI





Model selection

- Two-class classification example: AD vs HE





Display results & stat

- Example: AD vs HE classification




Compute weight maps




Display weights

MRI





PET



Cortical thickness



CADSuM
(Classification of AD Sufferers using Machine Learning)



- Surfstat toolbox has to be added to Matlab routines

Thickness view

- Note: the cortical thickness files has to be already calculated by free surfer!



Thickness stats

The data from a single subject or a group of subject can be analyzed.

An example of longitudinal analysis: multiple scans of one AD subject in 4 intervals.


RID	SITEID	VISCODE	DXCURREN
5	1	bl	2
5	1	m18	1
5	1	m36	1
5	1	m54	1

RID =
Subject ID



VISCODE =
time Interval

DXCURREN =
3 = AD
2 = MCI
1 = Healthy

Results indicates an improvement: Increased cortical thickness from interval 3 (m36) to interval (m54) e.g. at Parahippocampus area~ transition from MCI to Healthy conditions.





Feature analysis (stat toolbox)



Intracranial brain volumes (MRI)

Average intensity (PET)



Feature selection (statistical analysis)

- Statistical tests, t-test (2-group analysis) and ANOVA (more than 2 groups) are utilized to analyses features.
- This is under reconstruction!

References

1. Yousofzadeh V, McGuinness B, Maguire L and Wong-Lin K, 2016. A kernel-based Gaussian process for combined MRI and PET in Alzheimer's disease. In: Alzheimer's Research UK (ARUK) 2016 Conference, Manchester, UK. 8-9 March 2016.
2. Yousofzadeh V, McGuinness B, Maguire L and Wong-Lin K. Multi-kernel learning with Dartel enhances MRI-PET classification and prediction of Alzheimer's disease: group and individual analyses, Submitted to Front. Hum. Neurosci.